protect.Maternity belt

SANOMED

Anziehenleitung • Fitting Instructions • Mode d'emploi • Indicaciones para la colocación • Instruções de utilização • Brugsanvisning • Råd om användning • Bevestigingsinstructies • Istruzioni per indossare l'ortesi • Návod k nasazení • Οδηγίες εφαρμογής • Upute za oblačenje • Instrukcje zakładania • Рекомендации по надеванию • Kullanım talimatı

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2

protect.Maternity belt

Support belt for reducing back pain during pregnancy

Indications

Lower back pain during pregnancy

Contraindications/side effects

Currently none known

Mode of action

The integrated stabilisation rods in the lumbar region and the adjustable tension of the tension straps reduce the frontal strain that is caused by the increasing weight of the abdomen, relieving pressure on the abdominal wall, the pelvis and the spine.

Properties

- Velcro fastening, allows simple application, adjustment and removal
- Stabilisation rods in the lumbar region prevent the material from rolling
- Reduction of lower back and back pain through noticeable relief for back muscles
- Anatomical form ensures perfect fit.
- Proven support function

Instructions for putting on

- Place the maternity belt around your hips and position it so that the stabilisation rods are located over your spine in the middle of your lower back. Next fasten the two ends at the front under your abdomen (Fig. 1). Make sure that the belt is positioned well and does not bother you when walking or sitting.
- To fasten the elastic tension straps, place your hands in the loops of the two tension straps on the outermost side (Fig. 2). Pull these forwards until you reach the desired level of tension. Using the Velcro fastening, gently press the straps to attach them to the belt.

Washing instructions/care information

To maintain the adhesion of the tension straps and the fastening system and to prevent other damage, we recommend closing the fastening before washing and attaching the tension elements to the fastening components. Use a wool or fine wash detergent and wash at 30°, preferably as handwash. Avoid rubbing and spinning, rinse the material thoroughly and lay out flat to dry. Do not wring, spin, tumble-dry, dry on radiators or place in direct sunlight. Do not iron. Do not remove stains with benzene or other chemical cleaning agents. Do not dry-clean.

Storage

Store in a dry environment at room temperature. Protect against heat, direct sunlight and humidity.

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©medi Orthopädie