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

## The main characteristics of the elastic rope and how to attach it

Before beginning to exercise, attach the elastic rope and slip the grip over the hand.

Before exercising attach the elastic rope to a firm construction (for example, a heating pipe, a door handle or a table leg). At the other end,


The elastic rope is stretched only to half its maximum stretched length. Exercise is slow, using limited, easily managed force. The rope can be set to four different levels of force. If the rope is attached with the black extension, resistance of 1 kp is obtained. This force is the same as the force exerted on the body

## The main characteristics of the elastic rope and how to attach it

Before beginning to exercise, attach the elastic rope and slip the grip over the hand.

Before exercising attach the elastic rope to a firm construction (for example, a heating pipe, a door handle or a table leg). At the other end,

The elastic rope is stretched only to half its maximum stretched length. Exercise is slow, using limited, easily managed force. The rope can be set to four different levels of force. If the rope is attached with the black extension, resistance of 1 kp is obtained. This force is the same as the force exerted on the body

Producer:

Smíšek Richard SPS-EU s.r.o.  
Na Úbočí 2250/10, Prague 8  
Česká republika

Tel: 00420-284 810 231

E-mail: [spirstab@spirstab.com](mailto:spirstab@spirstab.com)

[www.spiralstabilization.com](http://www.spiralstabilization.com)

**IMPORTANT WARNING:**

We stretch the elastic rope only into half of its rest length.

In case of any sign of damage to its part, above all, it is about extension, it is necessary to stop exercising with a rope and this replace the damaged part.

Individual rope components can be ordered.

We will be happy to send them to you.

The rope is for exercise only and is necessary especially to prevent children from straining them and fired against other people's eyes

Producer:

Smíšek Richard SPS-EU s.r.o.  
Na Úbočí 2250/10, Prague 8  
Česká republika

Tel: 00420-284 810 231

E-mail: [spirstab@spirstab.com](mailto:spirstab@spirstab.com)

[www.spiralstabilization.com](http://www.spiralstabilization.com)

**IMPORTANT WARNING:**

We stretch the elastic rope only into half of its rest length.

In case of any sign of damage to its part, above all, it is about extension, it is necessary to stop exercising with a rope and this replace the damaged part.

Individual rope components can be ordered.

We will be happy to send them to you.

The rope is for exercise only and is necessary especially to prevent children from straining them and fired against other people's eyes